

Product Information Bulletin



Apparent viscosity, intrinsic viscosity, and molar mass, which are all indicative measurements of depolymerization, have been observed to change significantly upon time under certain storage conditions. Figure 2 shows example data for the apparent viscosity (measured in 1% (w/w) sodium alginate solution at 20°C using Brookfield viscosimeter at 20 rpm) from a stability study conducted on PRONOVA UP LVG (low viscosity, high-G sodium alginate). No significant change in apparent viscosity is apparent after 5 years of storage at 2-8°C/ambient relative humidity (RH), or at -18°C. For storage at 25°C/60% RH, apparent viscosity displayed approximately 20% decrease after 5 years, while at 40°C/75% RH, apparent viscosity was reduced by approx. 40% after 12 months.




Figure 2: Apparent viscosity values observed during 60 months of storage of PRONOVA UP LVG.

Based upon stability studies, NovaMatrix recommends refrigerated storage (2-8°C) of all PRONOVA UP sodium alginates. We label PRONOVA UP LVG/LVM/MVG/MVM with a shelf life of 3 years under refrigerated storage. PRONOVA UP VLVG and VLVM are new products currently not labeled with a shelf life, as stability studies are pending. We expect these very low viscosity products to be at least as stable as PRONOVA UP LVG/LVM/MVG/MVM. All PRONOVA UP sodium alginates are delivered in a closed polyethylene container. During storage, the lid of the container should be kept closed.

Our stability programs show that PRONOVA UP sodium alginates can be stored at room temperature (25°C/ 60% RH) for 6-9 months without significant reduction of apparent viscosity. Therefore, we do not consider cool/cold shipment of products to customers a necessity, as transportation and transit times are normally a matter of a few days.

PRONOVA UP specifications with respect to apparent viscosity are intended to be broad and suitable for a range of applications. We can offer customer-dedicated specifications if required.

References:

- Smidsrød, O., Haug, A., and Larsen, B. "Degradation of alginates in the presence of reducing compounds", *Acta. Chem. Scand.* **17** (1963), p2628-2637.
- Holme, H.K., Lindmo, K., Kristiansen, A., and Smidsrød, O. "Thermal depolymerization of alginate in the solid state", *Carbohydr. Polym.* **54** (2003), p431-438.
- Holme, H.K., Davidsen, L., Kristiansen, A. and Smidsrød, O. "Kinetics and mechanisms of depolymerization of alginate and chitosan in aqueous solution", *Carbohydr. Polym.*, Manuscript in press (2008).

FMC BioPolymer AS d/b/a NovaMatrix

Head Office, Operations, R&D, Order Center and Customer Service Norway:

Industriveien 33
N-1337 Sandvika
Norway

Tel.: +47 67815500
Fax: +47 67815510
E-Mail: novamatrix_info@fmc.com

Sales

Europe and Asia:

Tel.: +33 296541246
Cell: +33 613213039

Sales

Americas:

1735 Market Street
Philadelphia PA 19103
USA

Tel.: +1 215 299 6264
Fax: +1 215 299 6821

Technical Service

Europe and Asia:

Tel.: +47 67815500
Fax: +47 67815510
E-Mail: novamatrix_info@fmc.com

Technical Service

Americas:

Tel.: +1 609 963-6297
Fax: +1 609 963-6329
E-Mail: novamatrix_info@fmc.com

Patents

NovaMatrix/FMC Corporation does not warrant against infringements of patents of third parties by reason of any uses made of the product in combination with other material or in the operation of any process, and purchasers assume all risks of patent infringement by reason of any such use, combination, or operation.

Warranty

Because of the numerous factors affecting results, NovaMatrix/FMC products are sold under the understanding that purchasers will make their own tests to determine the suitability of these products for their particular purpose. The several uses suggested by NovaMatrix/FMC Corporation are presented only to assist our customers in exploring possible applications. All information and data presented are believed to be accurate and reliable, but are presented without the assumption of any liability by NovaMatrix/FMC Corporation.

Technical Service

The information contained in this bulletin is intended to be general in nature. Techniques and data pertaining to specific uses for NovaMatrix/FMC products and new developments will be published periodically in the form of supplemental application bulletins. Our technical staff is ready to offer assistance in the use of NovaMatrix/FMC products.

Regulatory Status

PRONOVA™ sodium alginate meets the standards set forth in the current editions of the United States Pharmacopeia/National Formulary and European Pharmacopoeia. PRONOVA™ sodium alginate satisfies ASTM F 2064 for use in tissue engineered medical products (TEMPs). PRONOVA™ sodium alginate is manufactured in compliance with cGMP and described in a DMF submitted to the US FDA.

Visit our web site at www.novamatrix.biz

NovaMatrix, FMC BioPolymer and PRONOVA are trademarks of FMC Corporation.

©2002 FMC Corporation
All rights reserved

